# **PRESS MASSAGE** slimming synergy **Revitalizing Action** PM **PRESS MASSAGE Draining Action PULSE MASSAGE Reducing Action PRESS WARM Tonifying Action** FM **MODELING PHYSIOSTIMUL**

OVERLINE



### LIGHTNESS TONE AND VITALITY

# slimming synergy












élite Resolution Warm & Pulse Massage s the innovative equipment by **Overline** for slimming the silhouette, using the synergy between heat and Press Massage, enhancing the benefits and aesthetic effects of both. The only system of its kind, a great novelty for the aesthetic sector. The effects of the **Overline** Press Massage are widely known. Its bodysuits and jacket with progressive inflation and overlapping internal air chambers, prevent the veno-lymphatic return, simulating the passage of the hands ensuring a uniform massage.

The **Overline** innovation is in combining the Press Massage to the slimming action of heat, a great combination that enhances the treatment, making it more effective and more comfortable. The final result of élite Resolution Warm & Pulse Massage is the progressive decrease of the circumferences, thanks to the reabsorption of the oedematous status and improvement of venouslymphatic circulation. The muds and bandages **iSol**, designed ad hoc for treatments with **Overline** equipment, have an important complementary action.





mm L480 x D440 x H1140 in L18,89 x D17,32 x H44,88



Memory Card 9465

Included 3130/12

DELUX 12 SECTORS
ARMS 4 SECTORS

Included 3130/8

**DELUX 8 SECTORS** 

Optional 3130/8WF

**ARMS 4 SECTORS** 



Press Massage contributes to the restoration of veno-lymphatic functions, improving the peripheral circulation. It has a draining effect, reducing swelling and muscle tension. It has an important detoxifying, revitalizing effect too, enhances the oxygenation of tissues and prevents the formation of cellulite \*. This treatment is very relaxing and generates a sense of well being and lightness.

### SUPER SILENT TECHNOLOGY

This apparatus is equipped with a new generation engine that makes the treatment silent and your work more efficient over time, with a life of 2500 hours, equal to 5000 treatments.

### 6 general body programs

(for 12 sectors equipment)

These treatments allow to work simultaneously on lower and upper limbs. Programs are diversified according to single needs and customer problems.

### 12 legs programs

8 sectors treatments that work on lower limbs. Programs are diversified according to single needs and customer problems.

### 4 arms programs

(optional armsuits for 8 sectors version)

Upward peristaltic treatment with 4 sectors that encourages the elimination of liquids, while reshaping the treated part.



The new PULSE system of élite Resolution Press Massage was created by the scientific laboratories Overline in order to make the drainage/massage treatment, performed by the equipment, more comfortable and more similar to a manual massage. The PULSE technology uses a computer system and managed by Resolution Press Massage, able to carry out a pulsed and continuous pressure in each sector.

The pulsed and intermittent pressure creates a high stress, improving the effects of the Press Massage treatment. The stimulating action of subcutaneous tissue improves the muscle contractility and metabolism, acting consequently on the arterial circulation and activating a natural drainage mechanism.



Heat generates pleasant and real beneficial effects on the body such as the increase of microcirculation with consequent improvement of interstitial exchanges and improvement of local metabolism. The increased metabolism has a trophic action, ie tissue nourishment, but also oxygenation, which acts effectively on imperfections caused by cellulite. The mechanism of exchange and activation of (local) metabolism is able to enhance detoxification. Finally, heat has important properties of activation of ingredients contained in Muds and specific Bandages iSol.



The shaping Physiostimulation uses frequencies ranging from 30 to 80 Hz. Ideal for an invigorating action recommended for those people with sagging tissues, this treatment acts on the entire length of the muscles, keeping skin elastic, stimulating the lymph flow in order to improve its transit, helping to release congestions created, further accelerating the outflow.

### Shaping treatment: preparation phases

The bands for Physiostimulation are applied over iSol bandages, with moisturizing, tonifying, anti-cellulite, slimming and draining effect.




### code 3130/12PWF



- Super Silent Technology
- 12 Sectors bodysuits
- Armsuits

### COMBINATIONS with 12 SECTORS bodysuits



### code 3130/8PWF



- Super Silent Technology
- 8 Sectors bodysuits
- Armsuits (Optional)

## COMBINATIONS with 8 SECTORS bodysuits



### Presso UPGRADE YOU ARE CREATING YOUR EQUIPMENT

Overline always thinks of increasing the work of professional aesthetic operators. Unfortunately, very often each of them has its own way of making aesthetics and would like "Her Technology": for this reason, thanks to the efforts of our laboratories, we have decided to offer the possibility of purchasing a basic equipment and then, based on individual needs, upgrade and improve it under different points of view and particularly based on the needs of your beauty center. You will no longer have the problem of having a device that you like but only half-way because you don't use all the functions, with Presso UPGRADE you can do whatever you prefer.







# TAKE CARE OF YOUR WORK DON'T WORRY ABOUT THE REST, OVERLINE IS THERE FOR YOU

**Overline** devices are guaranteed by defects in material and workmanship, ex-works, for a period of 36 months according to the General Conditions of Sale\*

The new Press Massage bodysuits are divided into feet, legs, undergarments with lymph node slits, and a side-adjustable bodice; sections concerning the arms start from the hand and are divided into 4 sectors, with a specific chamber for lymph nodes underarm. Bodysuits are adjustable in height and around the perimeter of legs and arms, with interchangeable inner air chambers.









# Enhance the benefits

# of Overline treatments



#### **FANGO-CELL**

Mud cream with white clay, greek hay and dead sea salts ISO.MUD.100 - 1000ml (Cabin)



#### **FANGO-SLIM**

Mud cream with green clay, active seaweed and green tea ISO.MUD.200 - 1000ml (Cabin)



#### LIPO-CELL GEL

Conductive gel for cellulite and adiposity with guaranà and green tea ISO.DREN.100 - 1000ml (Cabin)



### DRENA-TONE GEL

Conductive gel toning reactivating effect with gingseng ISO.TONE.200 - 1000ml (Cabin)



### CLEANSING SOFT MOUSSE

Dermocompatible cleanser Face & Body ISO.FOAM.150 - 150ml (Cabin & Home Care)



### JOJOBA FACE & BODY SCRUB

Gel scrub with Jojoba micro-granules ISO.SCRUB.200 - 250ml (Cabin)



#### HYDRABODY CELL

Body emulsion for cellulite ISO.BODY.100 Airless 200ml (Home Care)



### **HYDRABODY FIT**

Body emulsion toning effect ISO.BODY.200 Airless 200ml (Home Care)



#### DRAIN ACTIVE

Detox Draining Treatment ISO.BE.115

TONE UP Volume Effect Treatment ISO.BE.110

#### COMPACT FRESH

Tonifying Effect Treatment ISO.BE.105

### **REDU CELL**

Cellulite Reducing Effect Treatment ISO.BE.125

#### SLIM REDUCING

Slimming Effect Treatment ISO.BE.120

### THERMO CELL

Thermal Reducing Effect Treatment ISO.BE.100

**Overline** is a trademark **Technology srl**, a company that has obtained three important international certifications TUV RHEILAND. These certifications attest that the production chain complies with all regulations and international standards required for the production of biomedical equipment, electro-aesthetics and cosmetics, underlining how the excellence of Overline products is a competitive advantage for the beautician's work and for the satisfaction and tranquility of the end customer.





Occupational Health and Safety Assessment Series

OHSAS 18001:2007

www.tuv.com



Quality Management System

ISO 9001:2015

www.tuv.com ID 9105082302



Design, manufacture and service of electromedical devices

EN ISO 13485:2016

www.tuv.com D 9105082302

The manufacturer reserves the right to make modifications to the equipment at any time without prior warning. **Overline** equipment are designed and programmed to achieve results on beauty imperfections without exceeding the permitted aesthetic limits.

