OVERLINE

Precision LASER that Focuses on the Forms

Body Synergy Laser 980nm

Z IENNYHO



OVERLINE Body Synergy Laser 980nm


OVERLINE

19.**04**

BODY SYNERGY LASER 980nm

The brand new plates, with a slightly concave shape to better adapt to the body, allow the simultaneous and synergistic action of the best technologies for the remodeling of the body

ELECTRO REMODELING





1MHz Ultrasound - Electrolipolysis - LLLT

This combination proves to be one of the most effective and most powerful synergies to combat the localised adiposity and cellulite*. The Ultrasound acts selectively on the adipose cells and consequently on the fibrous component of the cellulite*, while the Electro-lipolysis acts on the adipose tissue, but especially on the lymphatic component and the microcirculation. The combination of these technologies with the LLLT makes it possible to act specifically on all the blemishes of cellulite* providing mutual reinforcement.

1MHz Ultrasound – Electro Remodeling – LLLT

The association between Ultrasound, Electro Remodeling and LLLT is a combination that generates an intense lipolysis action. Muscular contraction created by stimulation of the red muscle fibres makes it possible to collect the fats present in the tissues to tap into energy suitable for the contraction. This activity makes it possible to consume the lipids that otherwise are put into circulation and create new adipose deposits.

EFFECTIVE LIPOLYTIC AND REMODELING ACTION

The Low Level Laser Treatment (LLLT) technology includes a pre-treatment phase with electroporation that acts directly on the skin by creating electropores that, in relation to the length of the pulse, remain open from a few seconds to one or two minutes and permit a better permeability of the tissues in favour of a wide variety of hydrophilic molecules that otherwise would not be able to penetrate.

LLLT DREN

Programs that use a particular waveform with a "draining" effect thanks to the association of calibrated frequencies that permit, through a superficial muscle stimulus, the absorption of liquids.

The protocol provides for the use of **iSol Beauty** bandages

• ISO.BE.115 DRAIN ACTIVE (draining effect)

LLLT SLIM

In the Slim programs, stimulation is applied to the red muscle fibre assigned to energy consumption (aerobic) that improves the local metabolism, the "energy fuel" used to accomplish this job of training at the muscular level. The muscular stimulation consumes fats generating an improvement in the areas treated.

The protocol provides for the use of **iSol. Beauty** bandages

- ISO.BE.125 REDU CELL (reducing effect)
- ISO.BE.120 SLIM REDUCING (slimming effect)

LLLT LIPO

Programs that use the waveform and lipolytic frequency with continuous emission (without pauses) that acts at the level of the hypodermis, stimulating the division of the adipose cell into triglycerides > fatty acids > glycerol.

The protocol provides for the use of **iSol. Beauty** bandages

- ISO.BE.125 REDU CELL (reducing effect)
- ISO.BE.120 SLIM REDUCING (slimming effect)



LLLT TONIC

In the Tonic programs, the stimulation is applied to the intermediate muscle fibre that has mixed white fibre and red fibre characteristics. The stimulation of the local metabolism and the improvement in the vascularization and capillarization results in addition to an increase in muscle tone.

The protocol provides for the use of **iSol. Beauty** bandages ISO.BE.110 **TONE UP** (*volume effect*)

LLLT CELL EDEMA

The programs included in this context have a particular waveform that improves the purification and the drainage of excess liquids, creating a pleasant sensation of lightness.

The protocol provides for the use of **iSol. Beauty** bandages

- ISO.BE.125 **REDU CELL** (reducing effect)
- ISO.BE.115 DRAIN ACTIVE (draining effect)

*imperfection due to cellulite

The treatment at 360°

Thanks to this intense combination of energies, the treatment will offer to the Beautician very high performance in terms of result. Often when significant courses are performed on the imperfections of the localized adiposity and/or of the cellulite, the biggest problem occurs at the end of the courses where, frequently, one clashes with the atony is given by the sum of the treatments.

Overline has thought to provide synergistically to this unattractiveness, many times underestimated, but rather important, by inserting inside the equipment the Electro Remodeling technology. These absolutely pleasant stimulations at the perception and body level arouse the "muscle memory effect" creating filling and firming thanks to the toning of those areas where, because of unattractiveness, need to be reduced.

The shapes of the body will be redefined and the customers will immediately notice the changes in the terms of the remodeling of the silhouette, rediscovering the confidence of a physique which they were missing.





isol.Beauty Bandages, the perfect synergy

The blemishes caused by cellulite*, adipose deposits, lack of toning, orange peel skin, tissue pads and atony are the nightmare and the enemy of all women!

There are many methods to combat these problems. In addition to healthy eating and a correct life style, the synergy of the technologies offered by the Body Plate plates combined with the use of iSol Beauty bandages makes it possible to progressively and significantly expand the results.





Cuddles for Face

To regenerate, tone and elasticize the skin of the face, we have made available the handpiece for Face Radiofrequency so you can perform at the same time as the body treatments. A great advantage that translates into saving time and greater customer satisfaction, also optimizing the work for the beautician.

Body Synergy Laser 980nm

Brand New Overline technology for the remodelling and the synergy of different electrical waves to act in a targeted manner on the most frequent aesthetic imperfections that affect the body. Ideal for awakening the metabolism, improving the lymphatic circulation and microcirculation, draining, oxygenating, regenerating and stimulating the production of new Collagen, Elastin and Hyaluronic acid, deleting the extra rolls and the unattractive "orange peel" look, elite LLLT Resolution Body Synergy gives the skin firmness and elasticity.

The various programs that can be performed with elite LLLT Resolution Body Synergy exploit the use of innovative slightly concave plates, studied ad hoc to comply better with the shapes of the body, they act simultaneously thanks to a series of latest generation technologies: Low Level Laser (LLLT), Ultrasound, Electro-Lipolysis, Electro Remodelling and Electroporation, liquid retention and skin atony. The treatments generate a visible reduction of adipose tissue and the blemishes, skin firming. No recovery time necessary after treatment.

ref. 3245 Dimensions Inch. L21.65 x D18.90 x H45.30 Memory Card 9465





Enhance the benefits of Overline treatments

effect with Ginseng

(Cabin)



RF & STIM FACE EMULSION Face Emulsion for RF and Sti

Face Emulsion for RF and Stimol ISO.RF.100 - 250ml / 8,45oz. (*Cabin*)



DRENA-TONE GEL Conductive gel with toning reactivating

ISO.TONE.200 - 1000ml / 33,8oz.

(Cabin)

LIPO-CELL GEL Conductive gel for cellulite and adiposity with Guaranà and Green Tea ISO.DREN.100 - 1000ml / 33,8oz.

Contraction of the second seco

DRAIN ACTIVE Detox draining treatment ISO.BE.115

TONE UP Volume effect Treatment ISO.BE.110

*imperfections due to cellulite

REDU CELL Cellulite*reducing effect treatment ISO.BE.125

SLIM REDUCING Slimming effect treatment ISO.BE.120





TAKE CARE OF YOUR WORK, DON'T WORRY ABOUT THE REST, OVERLINE IS THERE FOR YOU

Overline devices are guaranteed by defects in material and workmanship, ex-works, for a period of 36 months (excluded the accessories subject to wear) according to the General Conditions of Sale*

Overline is a Technology S.r.l. Brand, a company that has obtained 3 important international TÜV Rheinland certifications. These certifications attest that the production chain complies with international standards and standards for the production of biomedical, electro-aesthetic and cosmetic equipment, underlining how the excellence of Overline products is a competitive advantage for the beautician's work and for satisfaction and the tranquility of the end customer





Occupational Health and Safety Assessment Series OHSAS 18001:2007





Quality Management System
ISO 9001:2015



Design, manufacture and service of electromedical devices

EN ISO 13485:2016

www.tuv.com ID 9105082302

Overline reserves the right to make changes to the equipment described, at any time and without notice. Overline equipment is designed and programmed to obtain results on imperfections without passing the allowed aesthetic limit.

www.tuv.com ID 9105082302

